

## Metabolic and Pharmacokinetic studies in peripheral tissues carried out by Microdialysis in chronically implanted conscious animals

### Why a new design of the Microdialysis Probe?

Microdialysis has become widely accepted as an *in vivo* technique for continuous monitoring of neurotransmitter release in both anaesthetized and conscious animals. However, applications of Microdialysis outside the CNS to study tissue metabolism, drug distribution and kinetics, have been less frequent until recently.


A number of papers where Microdialysis has been applied in different organs and body fluids was published. Most of these studies were done on anaesthetized animals, mostly because a Microdialysis Probe specially designed for body organs was not available. The **CMA/11** and **CMA/12 Microdialysis Probes** made for brain implantation are constructed of stainless-steel cannulae and hard plastic. They cannot be used

on conscious animals without causing pain, stress and additional tissue damage. Furthermore, a probe for peripheral tissues must have a special system for implantation, attachment to the tissue and connection to a perfusion pump.

These facts led to the construction of a new type of the Microdialysis Probe which is soft, flexible and easy to implant in virtually any organ in the body.

The **CMA/20 Microdialysis Probe** with a specially designed Introducer fulfils all the criteria for painless chronic use. The **CMA/20 Microdialysis Probe** together with the **CMA/120 System for Freely-Moving Animals**, provides a unique tool for continuous *in vivo* monitoring of molecules in a conscious animal for at least 24 hours.

### How is the CMA/20 Microdialysis Probe implanted?



*Fig. 1. The implantation procedure for the CMA/20 Microdialysis Probe (A) follows a similar procedure as used in clinical practice for intravenous catheterization. Briefly, a probe (A) is implanted into the tissue of an anaesthetized animal with help of an Introducer (B+C). The needle (C) is then removed and the CMA/20 Microdialysis Probe is inserted into the tubing (B). The plastic wing of the probe is sutured to the tissue. The tubing is removed by pulling and tearing it apart. The animal is placed in the CMA/120 System for Freely-Moving Animals and allowed to recover.*

## What are the benefits of using Microdialysis in peripheral tissues ?

The most fascinating feature of Microdialysis in conscious animals is that one can monitor

chemical changes in the body directly during a drug or food intake, stress, exercise, etc.




Fig. 2. The figure illustrates the very first application of the CMA/20 Microdialysis Probe to study the effects of acute stress (5 min tail pinch) on lactate, glycerol, glucose and urea levels in the fat of a rat. As seen, the 5 min stress (indicated by arrow) induced a dramatic increase in lactate, glycerol, and to certain extent glucose, within the first 30 minutes. Lactate returned quickly to basal levels, while glycerol and glucose remained increased for about another 60 minutes. Urea concentration is supposed not to be affected by activated

metabolism or blood flow and indeed, only a very slight increase in urea levels was measured. This fact may also serve as an evidence that recovery from the probe was not affected by the period of hyperactive movement of the rat during the stress. Thus it can be concluded that the observed changes of lactate, glycerol and glucose reflect the physiological response of the animal to stress, rather than experimental artifacts related to tissue movement.

## What are the requirements for a successful experiment?

The implantation of the Microdialysis Probe into tissues other than brain has some specific requirements. The most successful and reproducible experimental data are achieved using the CMA/20 Microdialysis Probes with 4 or 10 mm membranes and following recommended surgical

procedures. The CMA/120 System for Freely-Moving Animals should be used for connecting the animal to the CMA/100 Microinjection Pump and to the CMA/170, CMA/140 or the CMA/142 Fraction Collectors.




Fig. 3. The CMA/120 System for Freely-Moving Animals, CMA/100 Microinjection Pump, CMA/111 Syringe Selector and the CMA/170 Refrigerated Fraction Collector.

## What has been published about Microdialysis in peripheral organs ?

1. Gunnarson G., Jakobsson A.-K., Hamberger A., and Sjöstrand J. (1987) Free amino acids in pre-retinal vitreous space. Effect of high potassium and nipecotic acid. *Exp. Eye Res.* 44, 235-244.
2. Ben-Nun J., Cooper R.L., Cringle S.J., and Constable I.J. (1988) A new technique for in vivo intraocular pharmacokinetic measurements. *Arch. Ophthalmol.* 106, 254-259.
3. Ben-Nun J., Joyce D.A., Cooper R.L., Cringle S.J., and Constable I.J. (1989) Pharmacokinetics of intravitreal injection. Assessment of a gentamicin model by ocular dialysis. *Invest. Ophthalmol. Visual Sci.* 30, 1055-1061.
4. Nordenvall M., Ulmsten U., and Ungerstedt U. (1989) Influence of progesterone on the sodium and potassium concentrations of rat uterine fluid investigated by microdialysis. *Gynecol. Obstet. Invest.* 28, 73-77.
5. Hallström Å., Carlsson A., Hillered L., and Ungerstedt U. (1989) Simultaneous determination of lactate, pyruvate, and ascorbate in microdialysis samples from rat brain, blood, fat, and muscle using high-performance liquid chromatography. *J. Pharmacol. Meth.* 22, 113-124.
6. Dubey R.K., McAllister C.B., Inoue M., and Wilkinson G.R. (1989) Plasma binding and transport of diazepam across the blood-brain barrier. *J. Clin. Invest.* 84, 1155-1159.
7. Scott D.O., Sorensen L.R., and Iunte C.E. (1990) In vivo microdialysis sampling coupled to liquid chromatography for the study of acetaminophen metabolism. *J. Chromatogr.* 506, 461-469.
8. Caprioli R.M., and Lin S.-N. (1990) On-line analysis of penicillin blood levels in the live rat by combined microdialysis/fast-atom bombardment mass spectrometry. *Proc. Natl. Acad. Sci. USA* 87, 240-243.
9. Ståhle L., Segersvärd S., and Ungerstedt U. (1991) Drug distribution studies with microdialysis II: Caffeine and theophylline in blood, brain and other tissues in rats. in press.
10. Meirieu O., Pairet M., Sutra J.F., and Ruckebusch M. (1986) Local release of monoamines in the gastrointestinal tract: an in vivo study in rabbits. *Life Sci.* 38, 827-834.
11. Jarry H., Dietrich M., Barthel A., Giesler A., and Wuttke W. (1989) In vivo demonstration of a paracrine, inhibitory action of Met-enkephalin on adrenomedullary catecholamine release in the rat. *Endocrinology*, 125, 624-629.
12. Dorheim T.A., Wang T., Mentzer R.M.Jr., and Van Wylen D.G.L. (1990) Interstitial purine metabolites during regional myocardial ischemia. *J. Surg. Res.* 48, 491-497.
13. Van Wylen D.G.L., Willis J., Sodhi J., Weiss R.J., Lasley R.D., and Mentzer R.M.Jr. (1990) Cardiac microdialysis to estimate adenosine and coronary blood flow. *Am. J. Physiol.* 27, H1642-H1649.
14. Rosdahl H., Kehr J., Ungerstedt U., and Henriksson J. (1990) Muscle dialysis: Effects of locally administered adrenaline on extracellular levels of amino acids and lactate. *Acta Physiol. Scand.* 140, 1.
15. Fuchi T., Ungerstedt U., and Henriksson J. (1991) In vivo microdialysis of skeletal muscle. Measurement of changes in extracellular glucose concentration in the anaesthetized rat. submitted.
16. Lönnroth P., Jansson P.-A., Smith U., (1987) A microdialysis method allowing characterization of intercellular water space in humans. *Am. J. Physiol.* 253, E228-E231.
17. Arner A., Bolinder J., Eliasson A., Lundin A., and Ungerstedt U. (1988) Microdialysis of adipose tissue and blood for in vivo lipolysis studies. *Am. J. Physiol.* 255, E737-E742.
18. Bolinder J., Hagström E., Ungerstedt U., and Arner P. (1989) Microdialysis of subcutaneous adipose tissue in vivo for continuous glucose monitoring in man. *Scand. J. Clin. Lab. Invest.* 49, 465-474.
19. Lönnroth P., Jansson P.-A., Fredholm B.B., and Smith U. (1989) Microdialysis of intercellular adenosine concentration in subcutaneous tissue in humans. *Am. J. Physiol.* 256, E250-E255.
20. Hagström E., Arner P., Engfeldt P., Rössner S., and Bolinder J. (1990) In vivo subcutaneous adipose tissue glucose kinetics after glucose ingestion in obesity and fasting. *Scand. J. Clin. Lab. Invest.* 50, 129-136.
21. Arner P., Kriegholm E., Engfeldt P., and Bolinder J. (1990) Adrenergic regulation of lipolysis in situ at rest and during exercise. *J. Clin. Invest.* 85, 893-898.
22. Arner P., Kriegholm E., and Engfeldt P. (1990) In situ studies of catecholamine-induced lipolysis in human adipose tissue using microdialysis. *J. Pharmacol. Exp. Ther.* 254, 284-288.
23. Hagström E., Arner P., Ungerstedt U., and Bolinder J. (1990) Subcutaneous adipose tissue: a source of lactate production after glucose ingestion in humans. *Am. J. Physiol.* 258, E888-E893.
24. Jansson P.-A., Smith U., and Lönnroth P. (1990) Evidence for lactate production by human adipose tissue in vivo. *Diabetologia* 33, 253-256.
25. Jansson P.-A., Smith U., and Lönnroth P. (1990) Interstitial glycerol concentration measured by microdialysis in two subcutaneous regions in humans. *Am. J. Physiol.* 258, E918-E922.
26. Lönnroth P., and Smith U. (1990) Microdialysis - a novel technique for clinical investigations. *J. Int. Medicine* 227, 295-300.
27. Hargreaves K.M., and Costello A. (1990) Glucocorticoids suppress levels of immunoreactive bradykinin in inflamed tissue as evaluated by microdialysis probes. *Clin. Pharmacol. Ther.* 48, 168-178.
28. Kågedal B. Microdialysis of 5-S-cysteinyldopa in human melanoma tissue transplanted to athymic mice. Presented at 1st Int. Symp. on Microdialysis, 18-19 May, 1989, Indianapolis, U.S.A.
29. Anderson C., Andersson R., and Anderson T. Microdialysis a new bioanalytical sampling technique with relevant applications in skin contact reactions. Presented at 9th Int. Symp. on Contact Dermatitis, 17-19 May, 1990, Stockholm, Sweden.
30. Yokel R.A., Lidums V., McNamara P., and Ungerstedt U. (1991) Aluminium distribution into brain and liver of rats and rabbits following intravenous aluminium lactate or citrate: A microdialysis study. *Toxicol. Appl. Pharmacol.* 107, 153-163.

Note: The CMA Microdialysis Probes are not intended for use in humans. It is only suitable for laboratory research on animals.

**If you require further details on Microdialysis procedures, HPLC analysis, instrumentation or bibliography, please do not hesitate to contact:**

**CMA** *Microdialysis*

[www.microdialysis.com](http://www.microdialysis.com)

---

Box 2 • 171 18, Solna  
Tel 08-470 10 00 • Fax 08-470 10 50

73 Princeton Street, North Chelmsford • MA 01863 • USA  
Tel: +1 (978) 251-1940, (800) 440-4980, Fax: +1 (978) 251-1950