

- For early detection of brain ischemia
- "Gold Tip" visible on CT



# CMA 70

## Brain Microdialysis Catheter

This sterile, single use catheter is minimally invasive and designed for implantation in brain tissue. The dialysing membrane has been especially developed to achieve optimal diffusing characteristics. This allows a high recovery of substances from the extra cellular fluid into the catheter.

The membranes are available in 10 and 20 mm lengths, suitable for different target areas in the brain. The shaft is also available in different lengths making it possible to introduce the catheter by hand or stereotaxically. When introducing the catheter by hand the catheter is tunneled under the scalp with a tunneling needle. It is then inserted with the help of special forceps into the brain tissue through a hole drilled in the skull.

### The "Gold tip" makes the catheter visible on CT

The tip of the catheter contains a gold thread. The "Gold tip" is visible on CT-scanning and makes it possible to locate the exact position of the catheter.

### Early detection of local tissue ischemia


The CMA 70 Brain Microdialysis Catheter is connected to a CMA Microdialysis Pump. The outlet tube extends to a holder for a microvial that collects the dialysis sample. The microvials are usually changed every hour and the samples are transferred to and then analyzed in a CMA Microdialysis Analyzer. This provides an ideal clinical solution for early detection of ischemia in the brain.



**CMA** / *Microdialysis*

[www.microdialysis.com](http://www.microdialysis.com)

## Parts of the CMA 70 Brain Microdialysis Catheter



- |   |   |
|---|---|
| 1. Microvial - accessory (polystyrene + santoprene) | 9. Shaft (polyurethane)                 |
| 2. Luer lock connection (polycarbonate)             | 10. Dialysis membrane (polyamide)       |
| 3. Inlet tube (polyurethane)                        | 11. Gold thread within the catheter tip |
| 4. Vial holder (polycarbonate)                      | 12. Tunneling needle - accessory        |
| 5. Outlet tube (polyurethane)                       | 13. Forceps - accessory                 |
| 6. Stopper (silicone)                               |   |
| 7. Liquid cross (polysulfone)                       |   |
| 8. Protection tube (polyethylene)                   |   |



The distal part of the catheter has a gold thread (3x0.13mm) within the catheter tip, which makes the catheter location in the tissue visible on CT

## technical information

	material	length mm			Ø mm
		P000049	P000080	P000050	
<b>shaft</b>	polyurethane	60	60	100	OD 0.9
<b>membrane</b>	polyamide	10	20	10	OD 0.6
<b>inlet tube</b>	polyurethane	600	600	600	OD 1.0
<b>outlet tube</b>	polyurethane	220	220	220	OD 1.0

**CE** CE marked according to the Medical Device Directive, 93/42/EEC

**STERILE** Sterilised by  $\beta$ -radiation

**!** Storage temperature: 4-25 °C

**②** Single use only

## ordering information

	Ref. No.
<b>CMA 70 Brain Microdialysis Catheter</b> 4/pkg	<b>P000049</b>
<b>CMA 70 Brain Microdialysis Catheter</b> 4/pkg	<b>P000080</b>
<b>CMA 70 Brain Microdialysis Catheter</b> 4/pkg	<b>P000050</b>

### Accessories

Tunneling needle	P000055
Forceps	P000056
Microvials 250/pkg	P000001
Microvial Rack 12/pkg	P000028
Microvials in rack, Sterile 12x4	P000154



**CMA** / *Microdialysis*

[www.microdialysis.com](http://www.microdialysis.com)

Box 2, SE-171 18 Solna, Sweden  
Tel: +46 8 470 10 00 Fax: +46 8 470 10 50  
E-mail: [CMA@microdialysis.se](mailto:CMA@microdialysis.se)