



- Perfusion fluid for microdialysis in the brain
- Easy and safe to use

Perfusion Fluid CNS

for Microdialysis

Perfusion Fluid CNS is an isotonic sterile perfusion fluid especially developed for Microdialysis in CNS tissue. It is conveniently packed in glass ampoules, 5 mL, ready to use. Perfusion Fluid CNS is approved as a medical device according to the CE and the FDA medical directive.

Features	Isotonic, sterile	Package	Glass ampoules, 10 x 5 mL
Contents	NaCl 147 mmol/L KCl 2.7mmol/L CaCl ₂ 1.2 mmol/L MgCl ₂ 0.85 mmol/L Total Chloride contents: 153.8 mmol/L	Ref. No.	Perfusion fluid CNSP000151



CE marked according to the Medical Device Directive, 93/42/EEC



Storage temperature: 4-25 °C



Single use only



Sterilised by β -radiation



CMA / *Microdialysis*

www.microdialysis.com

Box 2, SE-171 18 Solna, Sweden
Tel: +46 8 470 10 00 Fax: +46 8 470 10 50
E-mail: CMA@microdialysis.se

US Office: 73 Princeton Street, North Chelmsford, MA 01863, USA
Tel: (800) 440-4980, (978) 251-1940
Fax: (978) 251-1950, E-mail: cma.usa@microdialysis.com